


# **MICRONYL** High performance nylon 66

- Easy integrity testable in situ
- Repeatedly steamable in situ and in autoclave
- Thermowelded construction
- EC-listed materials for Food contact
- FDA-listed materials per 21 CFR
- Bio-Safety per USP-Plastics
- Low filter extractables even with solvents
- Validation Guide available on request



MICRONYL cartridge includes Nylon 66 membrane at controlled porosity and provides high efficiency in bacteria retention.

The membrane is pleated with support and drainage layers in polyester which give high endurance versus thermal sterilization and hydraulic pulsation stress.

The intrinsically water wettability of nylon and polyester allows easy integrity testability. Typical application are downstream softening demineralization systems and in all process where cold sterilization is needed. Manufacturing is completed in a controlled environment; each filter is integrity tested.



Filter media	nylon 66
Upstream supports	polyester
Downstream supports	polyester
Internal Core	polypropylene
External Cage	polypropylene
End caps / Adapters	polyester

# **FOOD-SAFETY**

MICRONYL filter element materials meet (EU) regulation 10/2011 and its amendments, regulations (EC) 1935/2004 and 1895/2005.

#### **BIO-SAFETY**

Filter media and components pass USP CLASS VI Biological Reactivity and Chemical-Physical tests for USP plastics. Specific for "PS" grade: the filter meets USP "Water for injection" requirements for particle release and the effluent is Non-Pyrogenic per USP Bacterial Endotoxins (< 0,25 EU/ml).

QUALITY STANDARDS
Produced under a certified Quality System to guarantee traceability of manufacturing records and integrity testing results.

# **OPERATING CONDITIONS**


- max. continuous temperature	80 °C
- max. cumulative time of steam sterilization	20 hours at 121 °C / 13 hours at 125 °C / 5 hours at 135 °C
- sanitization with hot water	80 °C max
- sanitization with chemicals	Can be sanitized by standard chemical agents
- max. differential pressure	5,0 bar at 25 °C—2,5 bar 80 °C—0,3 bar 135 °C
- recommended change out differential pressure	2,0 bar at 25 °C
- recommended rinse up volume	3 liters/cartridge 10"


BACTERIAL RETENTION					
Brevundimonas d. Serratia m. Leuconostoc o. Escherichia c. Aeromonas h.				Aeromonas h.	
RATING CODE DY	≥ 10 <sup>6</sup>	≥ 10 <sup>11</sup>	≥ 10 <sup>11</sup>	≥ 10 <sup>11</sup>	≥ 10 <sup>11</sup>
RATING CODE DT	≥ 10 <sup>3</sup>	≥ 10 <sup>7</sup>	≥ 10 <sup>9</sup>	≥ 10 <sup>7</sup>	≥ 10 <sup>9</sup>
RATING CODE DK	/	≥ 10 <sup>6</sup>	≥ 10 <sup>8</sup>	≥ 10 <sup>6</sup>	/

# **INTEGRITY TESTABLE**




## WATER FLOW RATE FOR 10" CARTRIDGE





# MICRONYL MLE ORDERING INFORMATION

<u>DY</u> -



CODE
DY
DT
DK

CODE	
DY	
DT	
DK	

<u>BQ</u> -

<u>SB</u>

CODE	GASKETS	
No code	Standard	Silicone
E	On request	EPDM
V	On request	VITON
F	On request	FEP

CODE	PACKING TYPE
SB	Single box

CODE	NOMI NAL LENGTH
1	10"
2	20"
3	30″
4	40"

CODE	PRODUCT GRADE
BQ	Biological Grade; tested and prefluxed. Quality Certification in the box.
НІА	Product certified Halal BY WHA-IT00172-001

Data contained in this bulletin are informative and subject to change without notice. User is responsible for determining whether the product is fit for particular purpose and suitable for User's method of application.



DS-MLE-750-EN-20