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## INSTALLATION AND USE MANUAL FOR FILTER CAPSULES

### 1. INTRODUCTION

This manual contains the information necessary for the safe and efficient use and maintenance of Qualikap filter capsules.

Before performing any activity, the user is responsible for carefully reading the entire contents of this manual and must adequately train the personnel involved.

These filters are designed and manufactured in compliance with Directive 2014/68/EU PED Art. 4 Par. 3 according to the sound engineering practice (SEP). For this reason, the filter capsules do NOT bear the CE marking.

The filters has been designed for the uses described in this manual, any use other than those specified is considered improper and therefore potentially dangerous for people, property, and the surrounding environment.

### 2. WARNINGS AND LIMITS OF USE

It is the user's responsibility to evaluate the suitability of the filter capsules for a specific use and their adaptability to the process procedures.

The maximum operating conditions, based on the nature of the fluid and its hazardous nature, are as follows:

LIMITS OF USE		
Acc. To PED 2014/68/UE	GROUP 1	GROUP 2
LIQUID	3 barg @ 40 °C	5 barg @ 40 °C
GAS (*)	3 barg @ 40 °C	

(\*) For use with gas BEA Technologies recommends to use only the QKP-P and QKP-T capsules models.

Other capsules models (QKP-C, QKP-F, QKP-N, and QKP-S) are not suitable for use with GAS.

It is also recommended not to exceed the sterilization temperature indicated on the marking.

### 3. INSTALLATION

Before installation please follow the below recommendation:

- Carefully check the code and the filtration degree of the filter capsule;
- Extract the capsule from the plastic bag cutting only one side of it; do not use the capsule to force the opening of the bag;
- Orient the capsule so that the arrow printed on the marking and the process flow are in the same direction; check carefully the positioning of the o ring and secure it with the clamps first from one side than on the other side;
- during installation minimize the contact with the capsule in order not to alter the cleaning of the capsule itself;
- use the proper precautions if the capsule has been sterilized by Autoclave;
- pipes must be adequately supported in order to avoid to put additional load to the capsule.

### 4. STERILIZING

The sterilizing of the capsule has to be done by Autoclave as follows:

- position the capsule with the flow indicator arrow pointing down, to minimize condensate formation;
- **not autoclave capsules in horizontal position!!**
- the high and low point bleed valves must remain open during the steaming cycle by Autoclave;
- the in /out connectors must be protected using proper porous material to allow the vapor to reach the inside of the capsule; alternatively use one bag for Autoclave;
- do not connect the capsule with metallic part (clamp or piping) during the sterilization process;
- do not use the capsule if distort during the sterilizing process;
- if the filter is manufactured using polyethersulfone (PES) membrane and an integrity test has to be performed, it is recommended to wet completely the membrane before sterilizing;
- do not perform in line sterilization.

### 5. HYDROPHILIC MEMBRANE WETTING PROCEDURE

The procedure to wet the membrane inside the capsule is as follows:

- close both high and low point bleed valves positioned on the body of the capsule and close the valve located downstream (it is always suggested to have an upstream and downstream valve installed were a capsule is in place);
- slowly introduce the water (PW/WFI) or the process fluid, if required by the SOP, inside the capsule, open the bleed valves and let the liquid flow till the air bubbles are gone;
- close the high and low point bleed valves and at the same time open the valve downstream and adjust the water flow rate at 200 L/h per capsules;
- flush for at least 5 minutes

### 6. HYDROPHOBIC MEMBRANE WETTING PROCEDURE

Follow this procedure to wet the hydrophobic membrane:

- put the capsule inside a clean container with the flow indicator arrow pointing down;
- open bleed valves and fill the capsule with isopropyl alcohol (IPA/ water with a ratio of 60:40 vol./vol.) and allow the liquid to enter the membrane and valves;
- close the high and low point bleed valves and repeat 4-5 times the filling phase.

### 7. LIQUIDS FILTRATION

Capsules must not be subject to hydraulic shock; the pressurization must always be conducted following the flow indicator arrow (present on the body of the capsule).

Liquid filtration procedure is as follows:

- **slowly open the upstream valve** and let the liquid flow inside the capsule;
- **slowly open the high and low point bleed valves** and let the liquid flow till the air bubbles are over;
- open the valve downstream the capsule filter so that the filtration process is started;
- close the high and low point bleed valves;
- to minimize the eventual interaction between the capsule and the process fluid please flush the capsule with the first 2 liters of water (PW/WFI) or process fluid.

The capsules filter must be changed when the pressure drop is no more compatible with the process, in any case it cannot exceed the data reported in the Bea Technologies data sheet.

### 8. GAS FILTRATION

Please refer to the paragraph 5 liquid filtration procedure except for the vent phase that is no more necessary.

It is always necessary to keep high attention to the operations that lead the membrane to get wet and that will cause the liquid to be entrapped in the membrane pores (sterilization, steam-out), leading to an initial increase of the pressure drop, up to when the liquid will not be completely wiped out. The above mentioned pressure drop might lead to a distortion in the capsule structure lowering the filtration performances.

It is user responsibility to adopt all the precautions to avoid any distortion.

## 9. INTEGRITY TEST

The cGMP requires to test the filter element before and after the filtration process in order to evaluate the filter integrity.

The diffusion test has to be performed only with wet membrane (see point 4/5).

The hydrophilic membrane might be wet with water or if required by the SOP, with process fluid. The Hydrophobic membranes have to be wet with isopropyl alcohol (IPA/water with a ratio of 60:40 vol./vol.) or if required by the SOP, with process fluid.

The test parameters are reported in the Validation Guides available on demand. Bea Technologies is able to provide, on demand, the integrity test procedures.

User must adopt all suitable actions to prevent the process fluid to be contaminated by the liquid used in the integrity test.

## 10. DECOMMISSIONING AND DISPOSAL

Before disconnecting the capsule from the line, ensure that the system is depressurized or that the capsule is isolated upstream and downstream.

Drain the fluid inside and flush with inert gas.

Final disposal of the capsule must be carried out in accordance with the regulations in force in the user's country, respecting health and safety requirements.

## 11. STORAGE

The filter capsules must be kept in their original packaging, in a clean and dry space, without exposing it to the direct light, at a room temperature that varies from 5°C up to 35°C. The goods do not support heavy loads.

Bea Technologies recommend to use the product within the date indicated on the marking or on the certificate.

Bea Technologies declines any responsibility on the product if the recommended shelf life period is exceeded.

## 12. DISCLAIMER

Before any operation, ensure that the personnel responsible for operating the plant are adequately trained and equipped with appropriate protective systems for the hazards associated with pressure vessels and fluid leaks. Always consult the safety data sheets for the products being handled for information on how to manage any leaks.

For any further information and clarification do not hesitate to contact Bea Technologies at [info@bea-italy.com](mailto:info@bea-italy.com)